
KSK SECTIONAL MATRIX RING

- Rings used to hold the matrix in right place on the teeth.
- Simple nail design works with a variety of teeth. (\checkmark)
- Available with your clamp forceps. ()*Compatibility with KSK clamp forceps is guaranteed.
- Spring force is just right and easy to adjust.



ers : Dr. Noboru Tak Director of Takahashi Dental Office Clinical Lecturer of Tokyo Medical and Dental University





Specification Color: Black Autoclavable Material : Stainless Steel Made In Japan



KSK Sectional Matrix Rings 2 pcs Product Code : 06000010

KSK Sectional Matrix Rings set with Forceps Product Code : 06000020



Specification and Appearance are subject to change without prior notice for product improvement.

Manufacturer : DENTECH CORPORATION



