



*Ver.2025 - English*

DENTAL INSTRUMENTS    MADE IN JAPAN ●



# CONTENTS

## ENDODONTICS P.4

KSK Rubber Dam Clamps	P.4	RTP Reamers	P.21
Dental Dams	P.15	Finger Ruler	P.22
Rubber Dam Accessories	P.16	X-Ray Mesh Gauges	P.22
Rubber Dam Set	P.19	Resin Post Methacrylics	P.22
M Peeso Reamers	P.21	Retro Filling Retractor	P.23
Peeso Reamers	P.21	MTA carrier	P.23

## PROSTHETICS P.24

Matrix Bands	P.24	Tooth Contact Gauges	P.26
Matrix Retainer	P.24	New Inlay & Crown Setters	P.26
KSK Sectional Matrix Ring	P.25	Articulating Paper Holders	P.26
Tooth Separators	P.25	New Fine Lock Trays	P.27

## IMPLANTOLOGY P.28

Twist Drills Black	P.28	Drill Extenders	P.34
Step Drills For Guide	P.29	Drill Extension SHORT	P.34
Twist Drills For Guide	P.29	Technika Bone Mill Set	P.35
Twist Drills	P.30	Depth Gauges	P.35
Trephine Burs	P.31	Technika Repair Tap	P.36
Bone Trephines	P.31	Screw Remover	P.36
Circuler Knives	P.31	Gingiva Gauge	P.37
Surgical Burs	P.32	Bur Stands	P.37
Technika TK Drill System	P.33		
Multi-Adapter Wrench	P.33		
Implant Drill Checker	P.34		

## LABORATORY P.38

Vulcanite Burs	P.38	Air Jet Diamonds	P.38
----------------	------	------------------	------

## MISCELLANEOUS P.39

Lip & Cheek Retractor	P.39	pinlight2	P.42
Lip Expanders	P.39	Excavators	P.43
Big Oral Mirrors Real-Fine	P.40	Wax Spatulas 7A	P.43
Cotton Pliers	P.41	Evans	P.44
Tongue Guards & Saliva Ejector	P.41	Hollenback #3	P.44
Parallel Reflector	P.41	Woodson #1	P.44

### Autoclavable



Up to 121°C for 15 min



Up to 132°C for 5 min



Non - autoclavable



Made-to-order items

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS



## KSK Rubber Dam Clamps



### Features

- KSK Rubber Dam Clamps are made of specially selected and tempered stainless steel.
- Ideal spring action and high resistance against the corrosion.
- Various forms are available, so you can choose clamps in accordance with the teeth.
- 4000 times open-close metal fatigue test goes for every production lot to assure long-term spring action.
- For stable forms, our craftsmen adjust bending every KSK clamps one by one manually.
- All beaks are sharpened for best fitting on the teeth.
- "MATT" and "BLACK" can avoid eye fatigue of dentists from reflecting operation light.

### Finishing



#### ■ GLOSSY

Glossy finishing.



#### ■ MATT

Sand blasted non glossy finishing.




#### ■ BLACK


Not painted, Oxide layer finishing which is safe for human body

### Winged or Wingless

#### ■ WINGED



#### ■ WINGLESS



### [Important Safety Instructions]

It is likely to transform and to break when force is excessively added.  
Please do not expand rubber dam clamps more than the necessity.

KSK Rubber Dam  
Clamps


# Molar Clamps

Product Code

- Glossy
- Matt
- Black

3


- 01012003
- 01014003
- 01019003



For small lower molars with flat beaks.

W3


- 01012903
- 01014903
- 01019903



For small lower molars with flat beaks. Similar to 3 but wingless.

4

- 01012004
- 01014004
- 01019004



For upper molars with festooned beaks.

5


- 01012005
- 01014005
- 01019005



For upper molars with festooned beaks.

7


- 01012007
- 01014007
- 01019007



For big lower molars with flat beaks.

7A


- 01012907
- 01014907
- 01019907



For big lower molars with flat beaks.

8


- 01012008
- 01014008
- 01019008



For big upper molars with festooned beaks. Same size as 7.

8A

- 01012908
- 01014908
- 01019908



For upper molars with deeply festooned beaks which can grab subgingival part.

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS




## Molar Clamps

Product Code		Glossy
		Matt
		Black

W8A


 01012998  
 01014998  
 01019998



Similar to 8A but wingless.

12A


 01012912  
 01014912  
 01019912



Especially for the third upper left and lower right molars with flat serrated beaks, paired with 13A.

13A


 01012913  
 01014913  
 01019913



Especially for the third lower left and upper right molars with flat serrated beaks, paired with 12A.

14

 01012014  
 01014014  
 01019014



For partially erupted molars with deeply festooned beaks.


14A


 01012914  
 01014914  
 01019914



For partially erupted molars with deeply festooned beaks, bigger than 14.

15


 01012015  
 01014015  
 01019015



For big lower molars with flat and longer beaks, a strong spring.

18

 01012018  
 01014018  
 01019018



For upper & lower molars with flat beak.

24


 01012024  
 01014024  
 01019024



For upper left and lower right buccal molars' subgingival cavities. The buccal side has wider fringe slightly festooned beaks. Paired with 25.

25

 01012025  
 01014025  
 01019025



For lower left and upper right buccal molars' subgingival cavities. The buccal side has wider fringe slightly festooned beaks. Paired with 24.

# Molar Clamps

Product Code

Blue	Glossy
Grey	Matt
Black	Black

**26**

Blue	01012026
Grey	01014026
Black	01019026



For upper molars in case of exposure of the tooth neck with slightly festooned beaks. Available for the second bicuspid.

**28**

Blue	01012028
Grey	01014028
Black	01019028



For upper and lower molars with flat beaks, wide view.

**30**


Blue	01012030
Grey	01014030
Black	01019030



For lower left molar subgingival with slightly festooned beaks. The projection on the buccal beak slips into the furcation. Paired with 31.

**31**


Blue	01012031
Grey	01014031
Black	01019031



For lower right molar subgingival with slightly festooned beaks. The projection on the buccal beak slips into the furcation. Paired with 30.

**42**

Blue	01012042
Grey	01014042
Black	01019042



For upper left molars with flat beaks, paired with 43.

**43**

Blue	01012043
Grey	01014043
Black	01019043



For upper right molars with flat beaks, paired with 42.

**51**


Blue	01012051
Grey	01014051
Black	01019051



For buccolabial molar furcations with slightly festooned beaks. Enables to press down the gingiva.

**56**


Blue	01012056
Grey	01014056
Black	01019056



For lower molars with with festooned beaks.

**56A**

Blue	01012956
Grey	01014956
Black	01019956



For big lower molars with very slightly festooned beaks. One size larger than 56.

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS






## Molar Clamps

Product Code |  Glossy  
 Matt  
 Black

U67


01012067  
 01014067  
 01019067



Best for upper second molars with festooned beaks, compact size.

138


01012138  
 01014138  
 01019138



For upper right and lower left third molars with serrated offset beaks, paired with 139.

139


01012139  
 01014139  
 01019139



For upper left and lower right third molars with serrated offset beaks, paired with 138.

201


01012201  
 01014201  
 01019201



For upper molars with slightly festooned beaks.

202


01012202  
 01014202  
 01019202



For big lower molars with flat beaks. A little bigger than 200.

203


01012203  
 01014203  
 01019203



For lower right molars and the first bicuspid with flat beaks, Paired with 204.

204

01012204  
 01014204  
 01019204



For lower left molars and the first bicuspid with flat beaks, Paired with 203.

205


01012205  
 01014205  
 01019205



For upper molars with very slightly festooned beaks.

G-1

50580002  
 50550010  
 50550030



For lower left molars with cheek and tongue pressure shoving wings. No need to wear dental dam but gauze or cotton rolls. Paired with G-2.

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS




## Molar Clamps

Product Code | ■ Glossy | ■ Matt | ■ Black

G-2


■ 50580003  
■ 50550011  
■ 50550031



For lower right molars with cheek and tongue pressure shoving wings. No need to wear dental dam but gauze or cotton rolls. Paired with G-1.

DISTAL

■ 50030000



Available for molars, both deciduous and permanent easy to reach distal part of molar.


KSK Rubber Dam Clamps

## Bicuspid Clamps

Product Code  
■ Glossy  
■ Matt  
■ Black

0


■ 01012701  
■ 01014701  
■ 01019701



For narrow bicuspid cervical zone with flat beaks. Available for deciduous anteriors as well.

00

■ 01012702  
■ 01014702  
■ 01019702



For small bicuspid with flat beaks. Available for deciduous anteriors as well.

1

■ 01012001  
■ 01014001  
■ 01019001



For bicuspid root with festooned beaks.

1A

■ 01012901  
■ 01014901  
■ 01019901



For cervical cavities of bicuspid with festooned beaks.

2

■ 01012002  
■ 01014002  
■ 01019002



For lower bicuspid with flat beaks.

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS




## Bicuspid Clamps

Product Code |  Glossy  
 Matt  
 Black

2A

01012902  
 01014902  
 01019902



For bigger upper and lower bicuspid with flat beaks. A little bigger than 2.

W2A

01012992  
 01014992  
 01019992



Similar to 2A but wingless.

22


01012022  
 01014022  
 01019022



For upper and lower bicuspid with flat beaks, similar to 207 but wingless.

27

01012027  
 01014027  
 01019027



For upper and lower bicuspid with festooned beaks, similar to 206 but wingless.

29


01012029  
 01014029  
 01019029



For upper and lower bicuspid with flat beaks, for wider view.

206


01012206  
 01014206  
 01019206



For upper and lower bicuspid with festooned beaks. Similar to 27 but winged.

207


01012207  
 01014207  
 01019207



For upper and lower bicuspids with flat beaks. Similar to 22 but winged.

208


01012208  
 01014208  
 01019208



For big bicuspid. Similar to 207 but the fringe is little bigger.

209

01012209  
 01014209  
 01019209



For lower bicuspid with flat beaks. Similar to 208 but a little smaller and stronger spring.

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS

## Labial Clamps

Product Code |  Glossy |  Matt |  Black

**Y10** 01014999

For Pre-molar, Deep festooned beaks grab subgingival part, Wingless of No.1.

KSK Rubber Dam Clamps

### Labial Clamps

Product Code  
 Glossy  
 Matt  
 Black

**6** 01012006  
01014006  
01019006

For anteriors with slightly offset beaks.

**9** 01012009  
01014009  
01019009

For labial cavities on anteriors with slightly offset beaks.

**44** 01012044  
01014044  
01019044

For anteriors cuspid and cavities with flat beaks.

**90N** 01012090  
01014090  
01019090

Offset beaks and stronger spring.

**210** 01012210  
01014210  
01019210

For anteriors labial and canines and cavities on bicuspid labial.

**212** 01012212  
01014212  
01019212

For gum compression on anteriors with offset beaks.

**212 SA** 01018212  
01018213  
01018214

Similar to 212 but a little stronger pressure to gingiva.

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS



## Deciduous Clamps

Product Code |  Glossy  
 Matt  
 Black

KSK Rubber Dam Clamps

### Deciduous Clamps

Product Code  
 Glossy  
 Matt  
 Black

**54**  01015001

For upper right the first deciduous molars. (B/54)

**55**  01015002

For upper right the second deciduous molars. (A/55)

**64**  01015003

For upper left the first deciduous molars. (I/64)

**65**  01015004

For upper left the second deciduous molars. (J/65)

**74**  01015005

For lower left the first deciduous molars. (L/74)

**75**  01015006

For lower left the second deciduous molars. (K/75)

**84**  01015007

For lower right the first deciduous molars. (S/84)

**85**  01015008

For lower right the second deciduous molars. (T/85)

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY


MISCELLANEOUS

# Deciduous Clamps

Product Code	<div></div>	Glossy
	<div></div>	Matt
	<div></div>	Black

P-1


- 50580000
- 50550000
- 50550020



For lower left the first deciduous molars. (L/74)  
Paired with P-2.

P-2

- 50580001
- 50550001
- 50550021



For lower right the first deciduous molars. (S/84)  
Paired with P-1.

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISSCELLANEOUS



# ENDODONTICS

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS

KSK Rubber Dam  
Clamps

## Tiger Clamps

Product Code  
■ Glossy  
■ Matt  
■ Black

1T


■ 01012501  
■ 01014501  
■ 01019501



For bicuspid root with festooned serrated beaks.

1AT


■ 01012502  
■ 01014502  
■ 01019502



For bicuspids, especially for subgingival cavities with deeply festooned serrated beaks.

2T


■ 01012503  
■ 01014503  
■ 01019503



For lower bicuspids with flat serrated beaks.

2AT


■ 01012504  
■ 01014504  
■ 01019504



For upper and lower bigger bicuspids with flat serrated beaks. A little bigger than 2T

9T


■ 01012505  
■ 01014505  
■ 01019505



Used for labial cavities on anteriors with offset serrated beaks.

14T


■ 01012506  
■ 01014506  
■ 01019506



For partially erupted molars with serrated deeply festooned beaks.

56T

■ 01012507  
■ 01014507  
■ 01019507



For lower molars with festooned beaks.

## Dental Dams

### Dental Dams (Powder-Free )

#### Features

- KSK dental dams have a good resistance and elastic property.
- Pure natural rubber.
- "Thin" is made of color free material then capable for monitoring the condition of an patient during isolation treatment.

#### [Specification]

Packaging : 36 sheets / pack    Size : 6"x 6" (152mm x 152mm)  
 Thickness : apx 0.20mm(Medium), apx 0.18mm(Thin)  
 Color : Green(Medium), Ivory(Thin)  
 Scent : Non


Product Code 01023060

Product KSK Dental Dam Medium (cut)6"x6"

Product Code 01023070

Product KSK Dental Dam Thin (cut)6"x6"

Medium



Thin



ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISSCELLANEOUS

### Dental Dam Non-Latex

#### Features

- Used when there is a concern regarding an individual's possible allergic reaction to the proteins present in natural rubber.

#### [Specification]

Packaging : 20 sheets / pack    Size : 6" x 6" (152mm x 152mm)  
 Thickness : 0.20mm    Color : Blue    Scent : Non

Product Code 01023040

Product KSK Dental Dam (cut)6"x6" Non-Latex



### Dental Dam Non-Latex Powder-Free

#### Features

- Powder-free minimizes the risk of illness or injury caused by allergy for all dental health care personnel and patients.

#### [Specification]

Packaging : 15 sheets / pack    Size : 6" x 6"(152mm x 152mm)  
 Thickness : 0.25mm    Color : Purple    Scent : Mint

Product Code 01023050

Product KSK Dental Dam (cut)6" x6" Non-Latex Powder-Free







## Rubber Dam Accessories

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS

### KSK Rubber Dam Punch

[Specification]

Material : Stainless steel

Size : 153 x 70mm

Type : Ainsworth type.

Hole diameter : 1.4mm , 1.6mm , 1.8mm , 2.1mm , 2.3mm

Product Code 01021011

Product KSK Rubber Dam Punch



### KSK Rubber Dam Clamp Forceps

Features

Original design of DENTECH Rubber Dam Clamps are automatically off from the forceps at the place of operator fits them on tooth.

[Specification]

Material : Stainless steel

Size : 165 x 70mm

Product Code 01031040

Product KSK Rubber Dam Clamp Forceps



### Rubber Dam Frames Young Type

[Specification]

Material : Stainless steel

Adult : 108(w) x 95mm(h)


Pedo : 77(w) x 73mm(h)

Product Code 01022010


Product Rubber Dam Frame, Adult

Product Code 01022020

Product Rubber Dam Frame, Pedo



Adult



Pedo



### Rubber Dam Template

Features

For marking punch holes.


[Specification]

Material : Plastic

Size : 150 x 150mm

Product Code 08190020

Product Rubber Dam Template



## Clamp set

### Clamp set No.8

[Specification]

Material : Stainless steel Size : 90 x 90mm

With 8 clamps : #9 x 2, 2A, W2A, 5, 56, 203 & 204 or Clamp Organizer (Board only)

Product Code 01041000 Product Clamp Organizer

Product Code 01040000 Product Glossy

Product Code 01040080 Product Matt

Product Code 01040050 Product Black



Matt & Black



### Clamp set No.12

[Specification]

Material : Stainless steel Size : 150 x 90mm

With 12 clamps : #210, 9, 0, 2, 2A, 206, 207, 5, 56, 7, 201 & 202 or Clamp Organizer (Board only)

Product Code 01041010 Product Clamp Organizer

Product Code 01040010 Product Glossy

Product Code 01040090 Product Matt

Product Code 01040060 Product Black



### Clamp set No.18

[Specification]

Material : Stainless steel Size : 160 x 120mm

With 18 clamps : #210, 9, P-1, P-2, 22, 27, 29, 208, 209, 18, 138, 139, 7, 201, 202, 203, 204 & 205 or Clamp Organizer (Board only)

Product Code 01041020 Product Clamp Organizer

Product Code 01040020 Product Glossy

Product Code 01040100 Product Matt

Product Code 01040070 Product Black



ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS



## Deciduous Clamp set

[Specification]

Material : Stainless steel Size : 90 x 90mm

Complete kit of 8 patterns with special Clamp board.

#55-upper right E #54-upper right D #64-upper left D #65-upper left E  
#85-lower right E #84-lower right D #74-lower left D #75-lower left E

Product Code	01041061	Product	Clamp Organizer
--------------	----------	---------	-----------------

Product Code	01041070	Product	Matt
--------------	----------	---------	------



## Clamp set with cassette

[Specification]

Material : Stainless steel Size : 94 x 94 x 34 mm

Clamp set No.8 or Deciduous Clamp set with sterilization cassette

Product Code	01040001	Product	No.8, Glossy
--------------	----------	---------	--------------

Product Code	01040081	Product	No.8, Matt
--------------	----------	---------	------------

Product Code	01040051	Product	No.8, Black
--------------	----------	---------	-------------

Product Code	01041071	Product	Deciduous, Matt
--------------	----------	---------	-----------------



Rubber Dam Set

MTO

Matt & Black

132°C

121°C


Rubber Dam Technic set

[Specification]

Contents : ①Clamp set No.12  
②Rubber Dam Frame Large(Adult)  
③Rubber Dam Punch(Ainsworth)  
④Rubber Dam Clamp Forceps

Clamps : #210, #9, #0, #2, #2A, #206, #207, #5, #56, #7, #201, #202

Product Code	01030100	Product	Glossy
Product Code	01030120	Product	Matt
Product Code	01030110	Product	Black




Rubber Dam Complete set

[Specification]

Material : Stainless steel    Size : 274 x 184 x 34 mm

Rubber Dam Technic set with with sterilization cassette

Product Code	01030101	Product	Glossy
Product Code	01030121	Product	Matt
Product Code	01030111	Product	Black









## Rubber Dam Starter Kit

[Specification]

Basic assortment for rubberdam technique

Contents : ① Clamp set No.12,  
② Rubber Dam Frame(Adult)  
③ KSK Rubber Dam Punch(Ainsworth)  
④ KSK Rubber Dam Clamp Forceps  
⑤ KSK Dental Dam Non-Latex Blue(20 sheets)  
⑥ Rubber Dam Template

Product Code 01030130 Product Glossy



\*Dental Dam : not autoclavable

## Sterilization Cassette

Description

- Ultra light
- For washing, sterilization and storage

[Specification]

Material : Stainless steel

S : 94 x 94 x 34mm

L : 274 x184 x34mm

S for KSK Clamp set No.8 & KSK Deciduous Clamp set


L for KSK Rubber Dam Technic set

Product Code 50710605

Product Sterilization Cassette S

Product Code 50710610

Product Sterilization Cassette L



S



L

## M Peeso Reamers

[Specification]

Material : Stainless steel

Sharp shaped top

Pilot hole drills for root canal treatment.

Product Code 02312200 Product M Peeso Reamer CA #000 (set of 4 pcs)


Product Code 02312201 Product M Peeso Reamer CA #00 (set of 4 pcs)

Product Code 02312202 Product M Peeso Reamer CA #0 (set of 4 pcs)

Product Code 02312203 Product M Peeso Reamer CA #1 (set of 4 pcs)

Product Code 02312204 Product M Peeso Reamer CA #2 (set of 4 pcs)

Product Code 02312205 Product M Peeso Reamer CA #3 (set of 4 pcs)



## Peeso Reamers

[Specification]

Material : Stainless steel

Round shaped top

Pilot hole drills for root canal treatment.

Product Code 02312000 Product Peeso Reamer CA #000 (set of 4 pcs)


Product Code 02312001 Product Peeso Reamer CA #00 (set of 4 pcs)

Product Code 02312002 Product Peeso Reamer CA #0 (set of 4 pcs)

Product Code 02312003 Product Peeso Reamer CA #1 (set of 4 pcs)

Product Code 02312004 Product Peeso Reamer CA #2 (set of 4 pcs)

Product Code 02312005 Product Peeso Reamer CA #3 (set of 4 pcs)



## RTP Reamers

[Specification]

Material : Stainless steel


Round shaped

Final drills for root canal treatment.

Product Code 02312100  
Product RTP Reamer #1 (set of 4 pcs)

Product Code 02312101  
Product RTP Reamer #2 (set of 4 pcs)

Product Code 02312102  
Product RTP Reamer #3 (set of 4 pcs)



	L	φD	φd
# 1	9	1.35	0.92
# 2	10	1.65	1.17
# 3	10	1.95	1.48



## Finger Ruler

### Description


Can be attached on finger by free size. The adjustable stopper keeps it's position by spring.

[Specification]

Material : Stainless steel Measurement length : 32mm

Product Code 03050000

Product Finger Ruler



## X-Ray Mesh Gauge

### Description

Suitable for root canal treatment. Accurate examination before placing an implant. checking of calculus quantity and its area.

[Specification]


Cover Material : Polyestel Mesh Line Material : Stainless steel

●Dental Size

Size : 40 x 30mm Fine Grid : pitch 1mm Wide Grid : pitch 5mm

Product Code 07050010

Product X-Ray Mesh Gauge, Dental



## Resin Post Methacrylics

### Description

You use the core build upon by the indirect method. It disappears, becomes gas when casting. For the post of the same size and RTP reamer, post selection is easy.

[Specification]

Fine plastic burn - out posts.

Indirect only

Product Code 05000050

Product Resin Post Methacrylic #R1 (set of 50 pcs)

Product Code 05000051

Product Resin Post Methacrylic #R2 (set of 50 pcs)

Product Code 05000052

Product Resin Post Methacrylic #R3 (set of 50 pcs)

#R1(Blue)



#R2(Yellow)



#R3(Red)





# Retro Filling Retractor



Description

For reverse root canal treatment, Autoclavable

[Specification]

Material : Stainless steel    Measurement length : 176mm

Product Code	30600600
Product	Retro Filling Retractor #1



# MTA Carrier



Description

- Thin and long pipe suck up MTA and delivers into narrow and deep root canal.
- "Wide" suck up four times as amount as the existing product.
- Delicate Tip is designed to be replaceable in case of damaged.
- The handle is made of aluminum, light and making it easy to work with.

[Specification]

Overall length : 140mm

Tip outer diameter : Spare Parts φ0.8mm,Wide φ1.4mm

Material : Stainless (Carrier Parts), Aluminum (Handle)

Spare Parts contain Pipe and Extrude Bar 1each

Spare Parts Wide is optionally available.


Product Code	01060003
Product	MTA carrier

Product Code	01060013
Product	MTA carrier Spare Parts

Product Code	01060032
Product	MTA carrier Spare Parts Wide



Spare Parts Wide





# PROSTHETICS

## Matrix Bands




### Matrix Bands

[Specification]

Material : Stainless Steel Thickness : 0.04mm

Product Code	04021101	Product	Matrix Bands #A(set of 100 pcs)
Product Code	04021102	Product	Matrix Bands #B(set of 100 pcs)
Product Code	04021103	Product	Matrix Bands #C(set of 10 pcs)
Product Code	04021104	Product	Matrix Bands #D(set of 10 pcs)
Product Code	04021105	Product	Matrix Bands #E(set of 10 pcs)
Product Code	04021106	Product	Matrix Bands #F(set of 10 pcs)
Product Code	04021107	Product	Matrix Bands #G(set of 10 pcs)
Product Code	04021108	Product	Matrix Bands #H(set of 10 pcs)
Product Code	04021109	Product	Matrix Bands #I(set of 100 pcs)
Product Code	04021110	Product	Matrix Bands #J(set of 100 pcs)

※ Matrix Band #B-H is for sale while stock lasts.



### T Bands


#### Description

Can be used without the Matrix Retainer.


[Specification]

Material : Stainless Steel Thickness : 0.04mm

Product Code	04023111	Product	T Bands, Large (set of 100 pcs)
Product Code	04023112	Product	T Bands, Small (set of 100 pcs)



#### How to use T Bands



### Matrix Retainer #9-T Tofflemire set



#### Description

By a screw that has been improved, it will be a fine movement smooth.

[Specification]

Material : Stainless Steel Size : 64mm

Set : Matrix Retainers Tofflemire Type, Universal Adult, 9-T x 1pc, Matrix Bands #Ax 10pc

Product Code	04019300
Product	Matrix Retainers Tofflemire Type, Universal Adult, 9-T set



## KSK Sectional Matrix Ring



### Description

- Rings used to hold the matrix in right place on the teeth.
- Simple nail design works with a variety of teeth.
- Spring force is just right and easy to adjust.

[Specification]

Material : Stainless Steel

Product Code 06000010

Product KSK Sectional Matrix Rings 2 pcs

Product Code 06000020

Product KSK Sectional Matrix Rings set with Forceps



## Tooth Separators

### Elliot's Separator

### Description

The pushing screw can be attached both right and left side. Improved jaw design enable to be stable while pushing by the screw.

[Specification]

Material : Stainless Steel

Size : 37 x 67mm

Product Code 01071000

Product Elliot's Separator



### Ivory's Separator

### Description

With improved curved jaw design, thus separators become stable while pushing by the screw. It is applied in the front area of the upper and lower teeth.

[Specification]

Material : Stainless Steel

Size : 37 x 67mm

Product Code 01072000

Product Ivory's Separator





## Tooth Contact Gauges

### Description

For checking interdental contact points.  
In three colour coded thicknesses.

[Specification]

Material : Stainless Steel

Holder size:130 x 22 mm

Thickness: Blue 0.05 mm, Yellow 0.11 mm, Red 0.15 mm

Set: Holder x 1pc, Blue Gauges x 3pcs, Yellow Gauge x 1pc, Red Gauge x 1pc



Product Code 50050020

Product Tooth Contact Gauge Set

Product Code 50050021

Product Tooth Contact, Refill Gauge Blue (0.05mm)

Product Code 50050022

Product Tooth Contact, Refill Gauge Yellow (0.11mm)

Product Code 50050023

Product Tooth Contact, Refill Gauge Red (0.15mm)

## New Inlay & Crown Setters

[Specification]

Material : Stainless steel and silicone rubber

Long handle for molars has short pushing pin(SIZE : 79 x 24mm)

Short handle for bicuspid has long pushing pin (SIZE : 69 x 24mm)

Product Code 06092300

Product New Inlay & Crown Setter Set

Product Code 06092310

Product New Inlay & Crown Setter L Set

Product Code 06092320

Product New Inlay & Crown Setter S Set

Product Code 50050035

Product New Inlay & Crown Setter, Extra Rubber Pads 2 pcs



Silicone rubber :  
Dry heat are prohibited

※Replacement silicone rubber inserts 2pcs.

## Articulating Paper Holders

## Articulating Paper Holder, M-Type

### Description

It is made with milled in order not to move articulating paper and designed to be convenient for the operator.

[Specification]

Material : Stainless Steel

Size : 145mm

Product Code 06121000

Product Articulating Paper Holder, M-Type



( The inside )



# Articulating Paper Holder, Y-Type



Description


Sure holding of Articulating papers by pins and holes. Ideal to make for full dentures.

[Specification]

Material : Stainless Steel

Size : Large 74 x 60mm Small 66 x 48mm

Product Code	06122000
Product	Articulating Paper Holder, Y-Type, Large
Product Code	06122100
Product	Articulating Paper Holder, Y-Type, Small



# New Fine Lock Trays



Dry heat is prohibited

Description


"Retention rims" types. Brass , nickel plated . Complete with set of 8 size or individually each size .

Upper	U11	U12	U13	U14	mm
Wide	8 5	8 1	7 5	6 8	
Depth	5 9	5 7	5 6	5 4	
Height	1 9	2 0	2 0	2 1	

Lower	L11	L12	L13	L14	mm
Wide	8 0	7 6	7 6	7 2	
Depth	6 3	6 5	6 1	6 2	
Height	2 1	2 1	2 2	2 3	



Product Code	06040000
Product	New Fine Lock Tray Assort (set of 8 pcs)
Product Code	06041011
Product	New Fine Lock Tray #U11
Product Code	06041012
Product	New Fine Lock Tray #U12
Product Code	06041013
Product	New Fine Lock Tray #U13
Product Code	06041014
Product	New Fine Lock Tray #U14
Product Code	06042011
Product	New Fine Lock Tray #L11
Product Code	06042012
Product	New Fine Lock Tray #L12
Product Code	06042013
Product	New Fine Lock Tray #L13
Product Code	06042014
Product	New Fine Lock Tray #L14





## Twist Drills Black

### Description

A surgical drill with increased abrasion resistance, corrosion resistance, durability and hardness with DLC coating.

It comes with a movable stopper which allows for the length of the drill to be adjusted and a dedicated wrench. The adjustment range for the stopper is between 6 - 11mm for short, and 11 - 16mm for long.

### What's DLC?

DLC is short for "Diamond Like Carbon", and it is technology that coats a nano - level hard film with properties that resemble that of diamonds onto metallic surfaces. While having approximately 5 times the hardness compared to existing products, it has shown overwhelming low friction factor and it is recognized as the most technologically advanced among the many surface treatment techniques available.


[Specification]

Material : Stainless steel


Product Code	30073000	Product	Twist Drill "Black" Short 2.0S
Product Code	30073010	Product	Twist Drill "Black" Short 2.6S
Product Code	30073020	Product	Twist Drill "Black" Short 2.8S
Product Code	30073030	Product	Twist Drill "Black" Short 2.9S
Product Code	30073040	Product	Twist Drill "Black" Short 3.0S
Product Code	30073050	Product	Twist Drill "Black" Short 3.1S
Product Code	30073060	Product	Twist Drill "Black" Short 3.2S
Product Code	30073070	Product	Twist Drill "Black" Short 3.3S
Product Code	30073080	Product	Twist Drill "Black" Short 3.4S
Product Code	30073090	Product	Twist Drill "Black" Short 3.5S
Product Code	30073100	Product	Twist Drill "Black" Short 3.6S
Product Code	30073110	Product	Twist Drill "Black" Short 3.7S
Product Code	30073120	Product	Twist Drill "Black" Short 3.8S
Product Code	30073130	Product	Twist Drill "Black" Short 4.0S
Product Code	30073140	Product	Twist Drill "Black" Short 4.2S
Product Code	30073150	Product	Twist Drill "Black" Short 4.5S
Product Code	30073160	Product	Twist Drill "Black" Long 2.0L
Product Code	30073170	Product	Twist Drill "Black" Long 2.8L
Product Code	30073180	Product	Twist Drill "Black" Long 3.1L
Product Code	30073190	Product	Twist Drill "Black" Long 3.2L
Product Code	30073200	Product	Twist Drill "Black" Long 3.3L
Product Code	30073210	Product	Twist Drill "Black" Long 3.5L
Product Code	30073220	Product	Twist Drill "Black" Long 3.8L



Short



Long





# Step Drills For Guide




Description

Simplification of drill step has been made possible. It is ideal for pilot holes for tapered implants.

[Specification]

Material : Stainless steel

Product Code	30075000
Product	Step Drill "Black" for guide $\phi$ 2.4/2.8mm
Product Code	30075010
Product	Step Drill "Black" for guide $\phi$ 2.8/3.2mm
Product Code	30075020
Product	Step Drill "Black" for guide $\phi$ 3.2/3.6mm
Product Code	30075030
Product	Step Drill "Black" for guide $\phi$ 3.8/4.2mm
Product Code	30075040
Product	Step Drill "Black" for guide $\phi$ 4.2/4.6mm



# Twist Drills For Guide




Description

When using the 10mm surgical guide, it is marked to drill depths of 7mm, 8.5mm, 10mm.

[Specification]

Material : Stainless steel

Product Code	30074020
Product	Twist Drill "Black" for guide $\phi$ 2.0mm
Product Code	30074030
Product	Twist Drill "Black" for guide $\phi$ 3.0mm







# IMPLANTOLOGY

## Twist Drills




### Twist Drills External Cooling

#### Description

Our drills are made of stainless steel with high carbon. It assures enough hardness, good cutting efficiency and long life. But rather weak for rust. Careful maintenance will be required to avoid rust and to assure along life.

#### Good maintenance

- ① Immediately after use, drills should be dipped into anti-rust solution.
- ② Stopper should be removed while cleaning and autoclave.
- ③ Blood and bones should be removed completely after use.
- ④ Use only refined water for autoclave.
- ⑤ Use drill stand while autoclave, do not put into sterilization bag directly.



Product Code 30030200 Product Twist Drill Short 2.0S

Product Code 30030280 Product Twist Drill Short 2.8S

Product Code 30030310 Product Twist Drill Short 3.1S

Product Code 30030320 Product Twist Drill Short 3.2S

Product Code 30030330 Product Twist Drill Short 3.3S

Product Code 30030350 Product Twist Drill Short 3.5S

Product Code 30030380 Product Twist Drill Short 3.8S

Product Code 30030400 Product Twist Drill Short 4.0S

Product Code 30030420 Product Twist Drill Short 4.2S

Product Code 30030450 Product Twist Drill Short 4.5S

Product Code 30030550 Product Twist Drill Short 5.5S

Product Code 30030201 Product Twist Drill Long 2.0L

Product Code 30030281 Product Twist Drill Long 2.8L

Product Code 30030311 Product Twist Drill Long 3.1L

Product Code 30030321 Product Twist Drill Long 3.2L

Product Code 30030331 Product Twist Drill Long 3.3L

Product Code 30030351 Product Twist Drill Long 3.5L

Product Code 30030381 Product Twist Drill Long 3.8L

### Twist Drills Internal Cooling


#### Description

The drills with internal cooling system should be used also with external cooling. Drill speed should be controlled from 600rpm(thick drills) up to 1000rpm(for fine drills).

Cutting efficiency is coming bad by each use and bad drill will cause heat damages against bone. Only good drills should be used for operation. For very hard bones, it is recommended to use surgical burs prior to implant drills. Surgical burs will help easy drilling and long life. The fix screw of drill stopper is so small and it can not bear a strong force. Drills can be used also without stopper.

[Specification]

Material : Stainless steel



Product Code 30031200 Product Twist Drill Short 2.0 SI

Product Code 30031280 Product Twist Drill Short 2.8 SI

Product Code 30031310 Product Twist Drill Short 3.1 SI

Product Code 30031320 Product Twist Drill Short 3.2 SI

Product Code 30031330 Product Twist Drill Short 3.3 SI

Product Code 30031350 Product Twist Drill Short 3.5 SI

Product Code 30031380 Product Twist Drill Short 3.8 SI

Product Code 30031400 Product Twist Drill Short 4.0 SI

Product Code 30031420 Product Twist Drill Short 4.2 SI

Product Code 30031450 Product Twist Drill Short 4.5 SI

Product Code 30031201 Product Twist Drill Long 2.0 LI

Product Code 30031281 Product Twist Drill Long 2.8 LI

Product Code 30031311 Product Twist Drill Long 3.1 LI

Product Code 30031321 Product Twist Drill Long 3.2 LI

Product Code 30031331 Product Twist Drill Long 3.3 LI

Product Code 30031351 Product Twist Drill Long 3.5 LI

Product Code 30031381 Product Twist Drill Long 3.8 LI

Product Code 30032000 Product Screw for Implant Drill Stopper

Product Code 30032001

Product Hexagonal Wrench for Implant Drill Stopper

## Trephine Burs

### Description

Trephine Burs do not have blade at center, thus it will be rather difficult to get smooth cutting for a start. Enough exercise is recommended before use. After Trephine Burs penetrate into bone, then smooth cutting is available. To avoid choking on bone inside Trephine Bur, pull of Trephine Burs up and down, so that bone will always remain at patient. Remained bone can be collected by surgical curettes. Enough cooling and possibly low speed is recommended to avoid heat damage against bone.

[Specification]


Material : Stainless steel

Product Code 30060240 Product Trephine Bur 2.4

Product Code 30060270 Product Trephine Bur 2.7

Product Code 30060300 Product Trephine Bur 3.0

Product Code 30060330 Product Trephine Bur 3.3



Product Code 30060380 Product Trephine Bur 3.8

Product Code 30060430 Product Trephine Bur 4.3

Product Code 30060500 Product Trephine Bur 5.0

Product Code 30060600 Product Trephine Bur 6.0

## Bone Trephines

### Description

Bone Trephine is used for collection of big volume of bone. Enough exercise is recommended before use.

[Specification]

Material : Stainless steel


Product Code 30070500 Product Bone Trephine 5.0

Product Code 30070600 Product Bone Trephine 6.0

Product Code 30070800 Product Bone Trephine 8.0

Product Code 30071000 Product Bone Trephine 10.0

Product Code 30071200 Product Bone Trephine 12.0



## Circular Knives

### Description

For cutting soft tissue. Use with very low speed or by finger. Complete with finger grip for each piece.

[Specification]

Material : Stainless steel

Product Code 50010000 Product Circular Knife 3.3 without finger driver

Product Code 50010010 Product Circular Knife 3.5 without finger driver

Product Code 50010020 Product Circular Knife 4.0 without finger driver

Product Code 50010030 Product Circular Knife 5.0 without finger driver

Product Code 30051001 Product Finger Driver for Circular Knife





## Surgical Burs




### Surgical Bur External Cooling System

#### Description


- Round Burs  
Round burs for bone cutting, ruff cutting edges to prevent clogging. Full group blades ensures excellent cutting efficiency.

- Dr. Lindemann Bur  
Dr. Lindemann Bur has fully around blades.

Product Code	30011400	Product	Surgical Bur Round 141RF/010
Product Code	30011401	Product	Surgical Bur Round 141RF/014
Product Code	30011402	Product	Surgical Bur Round 141RF/018
Product Code	30011411	Product	Surgical Bur Round 141RF/025
Product Code	30011412	Product	Surgical Bur Round 141RF/031
Product Code	30011414	Product	Surgical Bur Round 141RF/040
Product Code	30011415	Product	Surgical Bur Round 141RF/050
Product Code	30011631	Product	Dr. Lindemann Bur 163RF/014
Product Code	30011641	Product	Dr. Lindemann Bur 164RF/018
Product Code	30011661	Product	Dr. Lindemann Bur 166RF/021



#### ● Dr. Lindemann Bur




### Surgical Bur Internal Cooling System

#### Description


These hully group surgical burs with internal cooling system enable very smooth cut and prevent heat damages for bone. Cooling fluid flow out blood thus operator can get clear view of operating spot. Prepare #40 K Reamers to clean internal cooling hold.

Product Code	30014711	Product	Surgical Bur Internal Cooling 471RF/023
Product Code	30011671	Product	Surgical Bur Internal Cooling 167RF/023
Product Code	30014731	Product	Surgical Bur Internal Cooling 473RF/031
Product Code	30014732	Product	Surgical Bur Internal Cooling 473RF/040
Product Code	30014733	Product	Surgical Bur Internal Cooling 473RF/050

#### ● Round Burs



#### ● Dr. Lindemann Bur





# Technika TK Drill System

### Description

Innovative dental drill for implant treatment.  
Safe and Secure drill systems  
Completed with a minimum 4drills, easy steps  
Ultra-low speed rotation, no water injection  
Autologous bone guide / Autologous bone collection

[Specification]

**Contents :**

- ① Surgical Bur Round 141RF/018
- ② Twist Drill Short 2.05 (External Cooling)
- ③ Technika TK Countersink (3.0, 3.5, 4.0, 4.5, 5.0)
- ④ Technika TK Drill (  $\varnothing$  3.0,  $\varnothing$  3.5,  $\varnothing$  4.0,  $\varnothing$  4.5,  $\varnothing$  5.0 )
- ⑤ Technika Surgical Cassette      \*final drill is not included



## IMPLANTOLOGY

<b>Product Code</b>	<b>30120000</b>	<b>Product</b>	<b>Technika TK Drill System</b>
---------------------	-----------------	----------------	---------------------------------

① Surgical Bur Round 141RF/018  
Material : Stainless Steel  
Length : 34 mm,  
Size : φ1.8mm  
Quantity : 1pc



②  
Twist Drill Short 2.0S (External Cooling)  
Material : Stainless Steel  
Length : 33 mm,  
Size :  $\phi 2.0$   
Quantity : 1pc



⑤  
Technika Surgical Cassette  
Material : Resin and Silicone  
Size : 99.5 mm x 143 mm x 61 mm  
Quantity : 1pc



③  
Technika TK Countersink  
Material : Stainless Steel  
Length : 24 mm,  
Size : 3.0, 3.5, 4.0, 4.5, 5.0  
Quantity : 1pc each



④ Technika TK Drill  
Material : Stainless Steel  
Length : 36 mm,  
Size :  $\phi 3.0$ ,  $\phi 3.5$ ,  $\phi 4.0$ ,  $\phi 4.5$ ,  $\phi 5.0$   
Quantity : 1pc each



- included drills



①      ②      ③

## Multi-Adapter Wrench

### Description

free-swing provides quick implantation.  
Conforms to major implant systems.

[Specification]

**Material : Stainless steel.**

Product Code	30080450	Product	Multi-Adapter Wrench
--------------	----------	---------	----------------------

Product Code	30080490	Product	Replacement Driver
--------------	----------	---------	--------------------

Product Code	30080500	Product	Square Adapter
--------------	----------	---------	----------------

Product Code	30080550	Product	HEX Adapter
--------------	----------	---------	-------------



## Multi-Adapter Wrench



Replacement Driver



Square Adapter



HEX Adapter



IMPLANTOLOGY | 33



# IMPLANTOLOGY

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS

## Implant Drill Checker

### Description


Can be measured diameter of 2.0 - 5.9mm.  
Can be measured length of 6.0 - 18.0mm.  
With length scale of 5 - 20mm. Can be used for drill stand.

[Specification]

Material : Stainless steel

Product Code 30080012

Product Implant Drill Checker





## Drill Extenders

### Description

Two types : External cooling/Internal cooling  
Extension length : External cooling is 16.5mm,  
Internal cooling is 18 mm.

[Specification]

Material : Stainless steel



Product Code 30080001

Product Drill Extension External Cooling

Product Code 30080002

Product Drill Extension Internal Cooling

## Drill Extension SHORT External Cooling

### Description

Can extend the drill length slightly 5mm.  
Easy attaching without tools.

[Specification]

Material : Stainless steel Diameter : 7.7 mm

Product Code 30079991

Product Drill Extension SHORT External Cooling



# Technika Bone Mill Set



Up to 121°C for 20 min  
Dry heat is prohibited

Description

- Simple and compact
- Easy operating  
The angle mount enables easy and smooth operating.
- Easy milling  
The unique shaped cutter with DLC coating enable easy and fine milling.
- Easy assembling  
The simple structure saves time for operating and cleaning.



Product Code	30080015
Product	Technika Bone Mill
Product Code	30080016
Product	Replacement Cutter


# Depth Gauges



Description

Graduations on 8, 10, 13, 15, 18, 20.  
Graduations can be exposed on X-ray film.  
Diameter : 1.3, 1.5, 1.9, 2.7, 3.0, 3.2, 3.4, 3.9mm.

[Specification] \_\_\_\_\_  
Material : Stainless steel



Product Code	30072130	Product	Depth Gauge 1.3
Product Code	30072150	Product	Depth Gauge 1.5
Product Code	30072190	Product	Depth Gauge 1.9
Product Code	30072270	Product	Depth Gauge 2.7

Product Code	30072300	Product	Depth Gauge 3.0
Product Code	30072320	Product	Depth Gauge 3.2
Product Code	30072340	Product	Depth Gauge 3.4
Product Code	30072390	Product	Depth Gauge 3.9





## Technika Repair Tap

### Description

- For internal threads repair.
- Made of high corrosion resistant stainless steel compatible with autoclave.
- Long lasting DLC coating.
- Safe Handle

### [Specification]

**Material :** Stainless steel (Tap and Handle)

**Surface Treatment :** DLC coating (Tap)

**Thread Diameter x Pitch :** M1.6 x 0.35 / M2.0 x 0.4

**Made in Japan**

Product Code	30011110	Product	REPAIR TAP M1.6 with Handle
Product Code	30011130	Product	REPAIR TAP M2.0 with Handle
Product Code	30011210	Product	M1.6 Tap only
Product Code	30011230	Product	M2.0 Tap only
Product Code	30011300	Product	Handle for Repair Tap

### ● Rwpair Tap with Handle(Basic Set)



### ● Repair Tap



For users having "Handle"

### ● Handle for Repair Tap



Compatible to both Tap

## SCREW REMOVER

### Description

- For removal of broken abutment screw.
- Made of high corrosion resistant stainless steel compatible with autoclave.
- Safe Handle

### [Specification]

**Material :** Stainless steel (Remover and Handle)

**Length :** S : 24mm / L : 32mm (For M1.6 and M2.0)

**Thread Diameter x Pitch :** M1.6 x  $\phi$ 1.1mm / M2.0 x  $\phi$ 1.4mm

Product Code	30010910	Product	M1.6 L with Handle
Product Code	30010915	Product	M1.6 S with Handle
Product Code	30010930	Product	M2.0 L with Handle
Product Code	30010935	Product	M2.0 S with Handle
Product Code	30011010	Product	M1.6 L Screw Remover only
Product Code	30011015	Product	M1.6 S Screw Remover only
Product Code	30011030	Product	M2.0 L Screw Remover only
Product Code	30011035	Product	M2.0 S Screw Remover only
Product Code	30011300	Product	Handle for Screw Remover

### ● Screw Remover with Handle(Basic set)



### ● Screw Remover



### ● Handle for Screw Remover



Compatible to both Tap




# Gingiva Gauge



**Description**

The first division shows 1mm when top of measuring block and needle are coincide. Wide of graduation is 1mm and pitch is 2mm, thus each graduation shows the gingival thickness of 2, 3, 4, 5, 6, 7mm. Put the measuring block in right angle to gingiva, then prick gingiva by light force until needle reaches to bone. Pull out carefully so that measuring block would not moves. Read the thickness at the back side of measuring block.

Product Code	30080013	Product	Gingiva Gauge
--------------	----------	---------	---------------



- Can be measured the thickness of gingiva from 2mm to 7mm in increments of 1mm during implant operation.
- For a more exact and detailed measure. It is better to use micrometer caliper separately.

# Bur Stands



**Description**

You can stand 6 burs.


[Specification]—

Material : Stainless Steel

Bur Stand 40(w70 x D40 x H15mm) Compliant bur length : Max 35mm

Bur Stand 45(W70 x D45 x H15mm) Compliant bur length : Max 40mm

Product Code	30080010	Product	Bur Stand 40
Product Code	30080011	Product	Bur Stand 45



ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS



## Vulcanite Burs

[Specification]

Material : Stainless Steel

Large and small types are available.

Large : Maximum diameter is 6 mm.

Small : Maximum diameter is 4.8 mm.



#AK



#BK




#CK



#FK



#A



#B



#C



#F



※ Straight hand piece type only.

Product Code	02101011	Product	Vulcanite Burs #A, LARGE
--------------	----------	---------	--------------------------

Product Code	02101012	Product	Vulcanite Burs #A, SMALL
--------------	----------	---------	--------------------------

Product Code	02101021	Product	Vulcanite Burs #B, LARGE
--------------	----------	---------	--------------------------

Product Code	02101022	Product	Vulcanite Burs #B, SMALL
--------------	----------	---------	--------------------------

Product Code	02101031	Product	Vulcanite Burs #C, LARGE
--------------	----------	---------	--------------------------

Product Code	02101032	Product	Vulcanite Burs #C, SMALL
--------------	----------	---------	--------------------------

Product Code	02101061	Product	Vulcanite Burs #F, LARGE
--------------	----------	---------	--------------------------

Product Code	02101062	Product	Vulcanite Burs #F, SMALL
--------------	----------	---------	--------------------------

Product Code	02111011	Product	Vulcanite Burs #AK, LARGE
--------------	----------	---------	---------------------------

Product Code	02111012	Product	Vulcanite Burs #AK, SMALL
--------------	----------	---------	---------------------------

Product Code	02111021	Product	Vulcanite Burs #BK, LARGE
--------------	----------	---------	---------------------------

Product Code	02111022	Product	Vulcanite Burs #BK, SMALL
--------------	----------	---------	---------------------------

Product Code	02111031	Product	Vulcanite Burs #CK, LARGE
--------------	----------	---------	---------------------------

Product Code	02111032	Product	Vulcanite Burs #CK, SMALL
--------------	----------	---------	---------------------------

Product Code	02111061	Product	Vulcanite Burs #FK, LARGE
--------------	----------	---------	---------------------------

Product Code	02111062	Product	Vulcanite Burs #FK, SMALL
--------------	----------	---------	---------------------------

## Air Jet Diamonds

### Features

Air Jet Diamond have air slits at cutting head, the slits supply strong air while cutting.

The strong air flow prevent the clogging of cutting materials. Thus, you can very smoothly and effective cutting.



Product Code	02091001
--------------	----------

Product	Air Jet Diamonds #1
---------	---------------------

Product Code	02091002
--------------	----------

Product	Air Jet Diamonds #2
---------	---------------------

Product Code	02091003
--------------	----------

Product	Air Jet Diamonds #3
---------	---------------------

Product Code	02091010
--------------	----------

Product	Air Jet Diamonds #1A
---------	----------------------

Product Code	02091011
--------------	----------

Product	Air Jet Diamonds #2A
---------	----------------------

#1A	
-----	--

Head Size

Φ1.2

#1	
----	--

Φ2.0

#2A	
-----	--

Φ3.1

#2	
----	--

Φ4.8

#3	
----	--

Φ6.0



## Lip & Cheek Retractor

### Features

To retractor cheeks and to have a better view.

### [Specification]

Material : Polycarbonate Set of 2pcs

Product Code 50570131 Product Lip & Cheek Retractor Set



Dry heat is prohibited

## Lip Expanders

### Features

Superb anterior access. Available in an adult, child or assorted(adult / child)set.

- Soft, plastic frame.
- Comfortable for the patient.
- Provides excellent anterior access.

### [Specification]

Material : Polypropylene

Product Code 50570010 Product Lip Expander Adult Size (set of 2 pcs)

Product Code 50570011 Product Lip Expander Child Size (set of 2 pcs)

Product Code 50570020 Product Lip Expander Adult & Child, (set of 2 pcs)



Dry heat is prohibited



ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS



# MISCELLANEOUS

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS

## Big Oral Mirrors Real-Fine

### Features


- Ideal for distortion - free intra - oral photos.
- Extensive rhodium coating on one side.
- Sturdy, durable.
- With detachable silicone handle.
- Autoclavable in the provided fabric bag.
- Reflection rate : 73%

#1 : Side views. #2 : Occulusal views (medium). #3 : Occulusal views (large).  
#4 : For lingual side views of posteriors.  
#5 : For lingual side views of anteriors. Handle for #4 , #5 : aluminum.  
Big Oral Mirror Set : Assort of #1,#2 and #3

Product Code	07000100
Product	Big Oral Mirror Set (#1~#3)
Product Code	07000101
Product	Big Oral Mirror #1
Product Code	07000102
Product	Big Oral Mirror #2
Product Code	07000103
Product	Big Oral Mirror #3
Product Code	07000104
Product	Silicone Handle for Big Oral Mirror
Product Code	07000105
Product	Big Oral Mirror #4
Product Code	07000106
Product	Big Oral Mirror #5
Product Code	07000109
Product	Aluminium Handle for Big Oral Mirror



Handle : Dry heat is prohibited



## Cotton Pliers

### Description


Available for picking up 0.01g.

### [Specification]

Coak type    Material : Stainless steel    Size : 162mm

Product Code    06301010

Product    Cotton Plier



## Tongue Guards & Saliva Ejector

### Description

Three different sizes metal reflectors are changeable and they protect tongue from turbine.

### [Specification]

Pipe : brass , chrome plated    Others : Stainless steel    Size : 138 x 80mm

Product Code    50050010

Product    Tongue Guards & Saliva Ejector



MOQ : 100pcs. Please contact us about the price every order.

## Parallel Reflector

### Features


Mirror with parallel lines.  
Autoclavable.

### [Specification]

Coak type    Material : Stainless steel  
Size : H 18 x W 48 mm, Thread Part M2.5 (P=0.45 L-5)

Product Code    50050100

Product    Parallel Reflector



MOQ : 100 pcs

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS



# MISCELLANEOUS

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS

## pinlight2

### Features

- The mirror and the light guide are detachable and autoclavable.
- The handle is IP X7 waterproof and capable of alcohol disinfection.
- Much brighter 90lm White LED Lighting used for detecting caries.
- Roder's HR FRONT mirror comes with size 4.

[Specification]

Light : 90LM, Cool white LED

Mirror : HR FRONT coating, size 4

Battery : Two AAA-size

Waterproof Performance IPX7 JIS C 0920 (IEC 60529)

Continuous lighting time : 90 min

Product Code	30080350
Product	pinlight2
Product Code	30080420
Product	Mirror holder, set of 4pcs.
Product Code	30080360
Product	Light guide
Product Code	30080370
Product	Replacement O ring, set of 3pcs.



Handle



Light guide  
Mirror top



Dry heat is prohibited  
Mirror holder



Mirror top No.4 included



# Excavators



Features

special hollow octagon handle which enables very light weight as well as ideal balance.  
high cutting performance and long life.

[Specification]  
Material : Stainless Steel Size : 165 mm x Handle : 5.8 mm

Product Code	50581100	Product	Excavators #19/20
Product Code	50581150	Product	Excavators #21/22

## #19/20



## #21/22



MOQ : 100 pcs

# Wax Spatulas #7A

[Specification]  
Material : Stainless Steel Size : 150 mm x Handle : 4.5 mm

Product Code	50581020	Product	Wax Spatulas #7A
--------------	----------	---------	------------------



MOQ : 100 pcs

ENDODONTICS

PROSTHETICS

IMPLANTOLOGY

LABORATORY

MISCELLANEOUS



# MISCELLANEOUS

ENDODONTICS

## Evans

[Specification]

Material : Spoon Stainless Steel, Handle Aluminum Size : 171 mm x Handle : 8 mm

Product Code 50581000

Product Evans



MOQ : 100 pcs

PROSTHETICS

## Hollenback #3

[Specification]

Material : Stainless Steel Size :166 mm x Handle : 5.8 mm

Product Code 50581040

Product Hollenback #3



MOQ : 100 pcs

IMPLANTOLOGY

## Plastic Filling Instruments, Woodson #1

[Specification]

Material : Stainless Steel Size : 174 mm x Handle : 5.8 mm

Product Code 50581060

Product Plastic Filling Instruments, Woodson #1



MOQ : 100 pcs

LABORATORY

MISCELLANEOUS

DENTECH Corporation (copyright, all rights reserved).

The Catalogue is intended to the practicing dentists, students of the medical universities and colleges, dental trade dealers. Design of the Catalogue, together with all the images, texts, graphic elements, text markers and identifiers are intellectual property of our company. For partial or full copying or loaning of any of the above mentioned elements for commercial purposes, you are obliged to ask for it and receive written permission.


Catalogue is for information only. The company does not carry any liability for possible miss prints, text mistakes or difference in products' colour or form of the items pictured on the images. The company leaves its' right for changes without prior consent.



# History of DENTECH

In 1907, Dentech Corporation, a pioneering manufacturer in Dental Instruments, was founded in Yanaka, Tokyo by Kenjiro Kimura, and it was originally named as Kimura Kousakusho.

We manufactured the first Japan made dental burs.



## After Second World War

In 1944, the first factory was destroyed due to Second World War. In 1947, after the war, we transferred our factory to Itabashi in Tokyo and reconstructed it.

However, in 1954, the strongest ever typhoon unfortunately struck around the area and the factory collapsed again.

In 1955, we rebuilt a new factory at the place nearby the destroyed factory to continue providing our products for the customers.



## First Advancement Overseas

After Second World War, a lot of dentists came to Japan from the U.S. accompanying the occupation army and needed trustworthy dental instruments with high quality.


Thus we were able to develop our products through the cooperation of them.

In 1971, we launched to export our products overseas and our first overseas customer was in the United States.

## Renamed as DENTECH

Kimura Kosakusho was renamed as Dentech Corporation in 1974. Since 1976 we have manufactured drills for implantology as the first implant drill manufacturer in Japan.

Nowadays most Japanese implant makers adopt our drills.



## The Current DENTECH

Today Dentech is well known as one of the companies which have the longest over 100 - year history in dental industry.

KSK rubber dam instruments, implantation drills and burs as well as many unique products with Japan-made high quality are well supported by dentists and professionals all over the world.



# Joint Performance with Japanese **CRAFTSMANSHIP** and Latest Machine Technologies




We possess CNC lathe, CNC Grinding machine, Heat treating machine, 3D CAD/CAM systems, and other various latest machines and equipments at our own factory in Japan. At the same time we are proud of the greatest technique by our craftsmen whose skill is unequalled with any others. We have an advantage in manufacturing extremely delicate dental/medical instruments using these two wheels which are identified by the machine automation and craftsmen's technique.

We have been developing, designing and manufacturing any kinds of dental and surgical instruments needed in clinical site around the world with outstanding quality for over a century.

It is a great honor for us that our products are highly respected all over the world. Our products would surely satisfy you.

Dentech has been assessed and certified as meeting the requirements of ISO 13485.

KSK Rubberdam system's products verified by the European medical device regulation, UK regulation and USA regulation.



*WEBSITE*



*e-CATALOGUE*